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
**SPITTOON DRAINING UNIT
GROUPE DE DRAINAGE DU CRACHOIR
SPUCKNAPF- DRÄNUNGSGRUPPE
UNIDAD DE DRENAJE ESCUPIDERA**

- MANUALE ISTRUZIONI
- OPERATOR'S HANDBOOK
- MANUEL D'UTILISATION
- GEBRAUCHSANWEISUNG
- MANUAL DE INSTRUCCIONES

CATTANI S.p.A.

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INDICE

Pagina

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INDEX

Page

— GENERAL RUNNING DATA.....	10
— INTRODUCTION	10
— SIGNS AND WARNINGS	10
— INSTALLATION	11
— FUNCTIONING	11
— RECOMMENDATIONS AND BIOLOGICAL HAZARD	11
— MAINTENANCE	11
— MAINTENANCE AND PERIODIC CHECKS BY AUTHORISED TECHNICIANS	12
— IMPORTANT NOTES	12
— TRANSPORT AND STORAGE	12
— TRANSPORT OF SECOND-HAND APPLIANCES	13
— ORDINARY MAINTENANCE RECORD	14
— EXTRAORDINARY MAINTENANCE RECORD	15

General running data

Spittoon Drainage Unit

Rated voltage:	24 V $\cong$ AC/DC $\pm 5\%$
Rated frequency:	50/60 Hz
Rated current:	0.2 A $\cong$ AC/DC
Absorbed power:	5 W
Insulation class:	I
Type of device:	B
Operating mode:	intermittent functioning
Liquid protection:	normal
Degree of protection against direct or indirect contact, type:	B 
Features of suction	
all-clear-relay contact:	Max 1 A - 24 V AC - Cos φ = 1
Flow rate	3 l/min

Introduction

The following presentation aims at illustrating the equipment and systems dealt with here in to professionals and repair experts. **We think it will also be useful to explain their operation and maintenance and inform** both experts and users about the hazards and precautions that should be taken to prevent accidents.

Signs and warnings

- **Biological danger, danger of infections from epidemic diseases.** 
- **High temperature.** 
- **General danger sign.** 
- **Compulsory direction of flow or rotation.** 

Signs cannot always fully express danger warnings, therefore it is necessary that the user reads the warnings and keeps them in due consideration.

Failure to respect a danger sign or warning may harm operator or damage the equipment. Do not remove protections, do not tamper with machines or their operation.

Installation (draw. page 40-41-42-43-44-45)

The spittoon drainage unit is to be **installed** inside the water system of the dental unit. It must **be placed vertically** and in a within-reach position allowing filter cleaning (fig. 2, part. A). Cylinder (part. C) can be assembled in the four positions, each one rotating by 90° with respect to filter holder (part. B).

The drainage tube of the spittoon shall be connected to the tube holder (part. D) of the filter holder (part. B). Discharge tube (part. E) connects the spittoon discharge unit to suction tube (fig. 1, part. F), located downstream the choking valve, by means of junction (part. G). **The unit operates at 24 V.** For electrical connections see electrical diagrams (fig. 4-5-6, page 43-44-45).

The little tube (part. N) must be connected to the hose (part. P) facing upwards. To ensure the proper functioning of the spittoon draining unit, if the filter (det. A) cannot be reached easily, the spittoon filter (det. R) must retain all solid particles up to 1.0 mm diameter.

Operation (draw. page 41)

For proper operation of the Spittoon Drainage Unit and of the Amalgam Separator, time about 20 seconds for rinsing of the Spittoon. The surgery staff shall be trained by the installer about the use and maintenance of the draining unit on the brand-new, non-contaminated appliance.

The water coming from the Spittoon enters the filter holder (part. B) and then the cylinder (part. C); when the water level reaches the highest probe (part. H), magnet (part. I) opens valve (part. L) and energises the relay, which starts the suction unit.

When the liquid level drops and the probe (part. M) does not sense liquid anymore, magnet (part. I) closes valve (part. L), while suction is timed to stop for 10 seconds with possible adjustments of up to one minute. When working on the machine, before disassembling the spittoon drainage unit start the aspirator and then pull upwards the cursor (detail Q) of magnet (detail I) to drain all the liquid inside the container (detail C).

Recommendations and biologic hazard

- **Biological danger, danger of infections from epidemic diseases.**

- **General danger sign.**

Before servicing any used equipment, carry out some washing operations with Puli-Jet plus (fig. page 12), or other disinfectant appropriate for this purpose. **Disconnect electric supply** and padlock mains switch if provided.

Before servicing the equipment wear disposable gloves, goggles and mask.

Aspirated debris is infected and contaminated from epidemic diseases, moreover the appliances to be serviced could be under pressure and the risk of contaminating splashes is clearly higher.

Maintenance (draw. page 41)

- **Biological danger, danger of infections from epidemic diseases.**


- **General danger sign.**

At the end of each working day, wash the system by sucking Pull-Jet diluted in warm water (50 °C) according to the instructions. Use Pulse Cleaner for the suction of the solution. This way the cleansing action of Puli-Jet plus will be accompanied by the necessary turbulence for a more thorough wash (see fig. below).

Pour the remaining part of Puli-Jet into spittoon. After cleaning as above mentioned, **take out the filter (fig. 2, part. A) while carefully avoiding any dripping; empty the contents of the filter as recommended for contaminated and toxic-noxious waste; then disinfect and rinse. If any solid particle goes beyond the filter (A), it is retained by the net (A1) placed at the base of filter housing (B). Aspirate these solid particles by means of the aspirator itself (see picture 1A page 41)**

Puli-Jet plus and Pulse Cleaner

Maintenance and periodic checks to be carried out by authorised technicians

- **Biological danger, danger of infections from epidemic diseases.**

- **General danger sign.**

Besides cleaning the filters, **carry out testing of magnet and electric circuits.** Bring rinsing of spittoon into operation and check the operation of magnet, the starting of suction with all terminals in the still position, and the status of relays and other electric contacts.

Important notes

- The manufacturer is at your disposal to supply spare parts, literature, instructions and anything that might be useful.
- Our concessionaires, agents, retailers, and authorised service engineers are provided with split up drawings, electrical diagrams, instructions and updating about maintenance and service.
- The guarantee and manufacturer's liability cease in case the appliances and/or the plants have been tampered with operations of any kind done by persons unfit and so unauthorised by the manufacturer.
- For any usage not contemplated or specified in this manual, please contact the manufacturer.

- On the web site: **www.cattani.it** you can find our **up-dated** instructions manuals. We recommend reading them, particularly for the **safety** updating.

Transport and storage

- Packed equipment can be transported and stored at a temperature range of - 10 °C + 60 °C.
- Packages must be kept away from water and splashing and cannot tolerate humidity >70%.
- Packages with the same weight can be stored in piles of three only.

Transport of second-hand appliances

- Prior to packing, clean and disinfect with Puli-Jet plus (see "Recommendations and biologic hazard" and "Maintenance").
- Place the appliance into a polyethylene bag, seal it and pack in a 3-layer corrugated board.

ORDINARY MAINTENANCE RECORD

Aspiration system

Record to be filled in at every ordinary maintenance operation

The ordinary maintenance must be entrusted to the surgery staff purposely trained. Before any maintenance operation, wear disposable gloves, mask, overalls, and goggles.

Date of maintenance Operator's signature		
Every evening: aspirate Puli-Jet plus with hot water by means of Pulse-Cleaner	YES	NO
Before getting near the appliance, cut the power off	YES	NO
Every day: replace filters on the dental unit and on the aspirator; put disinfecting antifoaming tablets into clean filters	YES	NO
Every three months: check the functioning of sound and visual signals	YES	NO
At intervals: check temperature in the technical room (from +5 to + 35 °C)	YES	NO
At intervals: check sound pressure level	YES	NO
Keep the technical room clear of anything that is not related to the machines contained therein, taking special attention of flammable material; make sure that there is no risk of formation of corrosive, flammable or explosive mixes	YES	NO
At every maintenance operation: check there is no hindrance (dust or other) to the proper ventilation of the machines	YES	NO

Photocopy this record and keep it after filling in

EXTRAORDINARY MAINTENANCE RECORD

Aspiration system

Record to be filled in at every extraordinary maintenance operation

Extraordinary maintenance is to be entrusted to a trained engineer provided with original spare parts. Before any maintenance operation, wear disposable gloves, mask, overalls, and goggles.

Photocopy this record and keep it after filling in

ENGLISH

Date of maintenanceOperator's signature		
Before getting near the appliance, cut the power off	YES	NO
At the first check up: check the clamping of terminals in the control panel	YES	NO
At the first check up: check the tension of the belt, if there is one	YES	NO
At the first check up: check the temperature of the technical room (from +5 to +35 °C)	YES	NO
Every three months: check the filters that the surgery staff can not reach	YES	NO
Every three months: check if there is any bacterial coating in the piping	YES	NO
Every three months: check the functioning of compensating valves and of thermal compensating ones	YES	NO
Every three months: check the functioning of electro-valves and of draining pumps	YES	NO
Every three months: check the functioning of probes and waste outlets	YES	NO
Every three months: check the tension of the belt, if there is one	YES	NO
Every three months: verify the sound and visual signals	YES	NO
Every six months: check the noise. The noise level (in dB (A) of each machine, which has been tested according to ISO 3647 (E) standard, can be found in instructions manuals and in our web site www.cattani.it	YES	NO
Every six months: check relays and remote switches contacts	YES	NO
Every six months: check the capacity of the capacitors	YES	NO
Every six months: test aspiration by checking flow and head; in case of reduced performances, take necessary measures.	YES	NO
A too loose belt is noisy at starting and stopping, and vibration can be seen in the loose section. A too taut belt makes the bearings whistle shrilly.		

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CONNECTION OF SPITTOON DRAINING UNIT TO ELECTRICAL CONTROL PANELS, IN SEMI-WET AND WET ASPIRATION SYSTEMS FITTED WITH ELECTRIC-PNEUMATIC CONTROL PANELS AND WITH ONE OR MORE ASPIRATION UNIT

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CONNECTION OF SPITTOON DRAINING UNIT TO ELECTRONIC CONTROL PANELS AND PERIPHERAL UNITS IN SEMI-WET AND WET ASPIRATION SYSTEMS



Fig. 5
Draw 5

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AC 33 SPITTOON DRAINING UNIT LEVEL CONTROL

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Fig. 6
Draw 6

ITALIAN PATENTS OR PATENT APPLICATIONS:

CATTANI: 1234828 - 1259318 - 1.187.187 - 1253460 - 233634 - 2337706 - 1294904 - 1305443 - 1310808 - 1340812
ESAM: 1225173 - 1253783 - 0791751

FOREIGN PATENTS OR PATENT APPLICATIONS:

CATTANI: US 4,787,846 - US 5,039,405 - US 5,002,486 - US 4,684,345 - US 5,330,641 - CH 0211808 - DE 0211808 - FR 0211808 - GB 0211808 - SE 0211808 - DE 0335061 - ES 0335061 - FR 0335061 - GB 0335061 - AT 0557251 - DE 0557251 - ES 0557251 - FR 0557251 - GB 0557251 - DE 0638295 - DK 0638295 - ES 0638295 - FR 0638295 - GB 0638295 - NL 0638295 - SE 0638295 - US 6,083,306 - US 6,090,286 - US 6,022,216 - US 6,423,124
ESAM: US 4,948,334 - DE 0351372 - ES 0351372 - FR 0351372 - GB 0351372 - US 5,779,443 - CH 0791751 - DE 0791751 - ES 0791751 - FR 0791751 - GB 0791751 - PT 0791751 - ES 107358 - FR 222.394/395

PENDING PATENT

CATTANI: EP 99830010.7 - EP 99830011.5 - EP 99830250.9 - EP 00830491.7 - EP 02425482.3
ESAM: EP 02425317.1 - MO 2001A000176 - US 2003-0044295

CATTANI s.p.A.

VIA NATTA, 6/A - 43100 PARMA - ITALY

TEL: +39 0521 607604 - SALE DEPT. FAX: +39 0521 607628

PURCHASING DEPT. FAX: +39 0521 607855 - ACCOUNTING DEPT. FAX: +39 0521 399966

<http://www.cattani.it> Email: cattani@tin.it

Company with Quality System Certified by DNV UNI EN ISO 9001/2000 - UNI CEI EN ISO 13485

esam s.p.A.

VIA NATTA, 4/A - 43100 PARMA - ITALY

TEL: +39 0521 607613 - SALE DEPT. FAX: +39 0521 607628

PURCHASING DEPT. FAX: +39 0521 607855 - ACCOUNTING DEPT. FAX: +39 0521 399966

<http://www.esam.it> Email: esamsa@tin.it

Company with Quality System Certified by DNV UNI EN ISO 9001/2000