REPROCESSING INSTRUCTIONS

"I FEEL
GOOD"

EMS+

EN: REPROCESSING INSTRUCTIONS

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BEFORE USE

Please read these Operating Instructions carefully as they explain all the most important details and procedures. Please pay special attention to the safety precautions. Always keep these Instructions close at hand.

 These Instructions are only applicable to the equipment they were delivered with.

To prevent injury to persons or damage to property, please observe the corresponding directives and symbols. This instruction follows the recommendations of RKI - Robert Koch Institut, Germany.

Translations in this manual are based on the English version. Thus the English text will be referred to in case of discrepancy.

DEVICES AND HANDPIECE CORDS

 The following instructions have been validated for all reusable products¹ of EMS dental devices. It remains the responsibility of the processor to ensure that the processing, as actually performed using the equipment, materials and personnel in the processing facility, achieves the desired result. This requires validation and/or validation and routine monitoring of the process. Likewise any deviation by the processor from the instructions provided should be properly evaluated for effectiveness and potential adverse effects.

AIRFLOW® and PIEZON® devices

Handy

PIEZON® Kit on your unit

CLEAN

Wipe the surface with a wet cloth until all visible contamination is removed.

DISINFECT

Wipe the surface with 70% isopropyl alcohol or a disinfection wipe with proven efficacy.

 Follow carefully the instructions provided by the disinfection solution manufacturer.

 **Do not use CaviWipes™.** It damages EMS products.

INTRODUCTION

EMS recommends compliance with cleaning, disinfection, packaging for sterilization and sterilization procedures according to ISO 17664.

 Always report adverse events related to device reprocessing directly to EMS.

 Products must be cleaned, disinfected and, if applicable, sterilized prior to first use. Do not reprocess the products over the authorized number of sterilization cycles. To replace, refer to the "Service life" section of these Instructions.

 Concentrations and contact times specified by the manufacturer of the cleaning and disinfection agent must be observed.

 Remember that sterilization cannot be achieved unless assembly components are cleaned and disinfected first.

If parts of these instructions are not clear or seem to be inadequate, do not hesitate to contact/inform EMS.

 **Users shall also observe any legal requirements applicable in their country as well as the hygiene regulations in effect at the hospital or clinic.**

This applies especially with regard to additional requirements for the inactivation of prions.

SERVICE LIFE

The products have been designed for a large number of sterilization cycles. The materials used in their manufacture were selected accordingly. However, with every renewed preparation for use, thermal and chemical stresses will result in ageing of the products.

 Always replace products that present signs of wear or early degradation, regardless of the number of sterilization cycles left unused.

 **DO NOT** expose the products to temperature exceeding 138°C.

 EMS handpieces must not be lubricated or otherwise exposed to oil whether directly or indirectly.

Lifetime : see product's Instructions for Use.

SYMBOLS

 Caution / Danger
Risk of property damage or injury

 Please note
Useful additional information and hints

 Prohibited

 Wear personal protective equipment

 Refer to instruction manual

 Sterilizable at up to 135°C in the autoclave


 Suitable for thermal disinfection

 Required time

Abbreviations used: min (minute) ; s (second)

¹ Hereinafter referred to as "products".

PROCESS

✓ PREPARE

✓ 20 s - Immediately after use, rinse the lumen(s) line of the handpiece/instrument with water. Coarse soiling must be removed immediately after application.

Safely transport to the reprocessing area to avoid any damage to the products, contamination of the environment or of the persons involved in the reprocessing process.

⚠️ ⌚ Cleaning is to be performed within 1 hour from use.

🧤 🧴 Wear personal protective equipment according to the type of preparation.

! AIRFLOW® and PERIOFLOW® Handpieces:

Always carry out handpiece powder unclogging and check for both lumens (water and powder) clearance before proceeding. Use Easy Clean.

! PIEZON® Endochuck:

Remove any installed file.
Remove the screw and separate the small O-ring.

! PIEZON® Handpieces:

Remove any installed instrument.
Remove the nose cap and separate the light guide and O-ring.

! PIEZON® Instruments with CombiTorque®:

Separate the PIEZON® Instrument and the CombiTorque® but keep each original instrument and CombiTorque® together throughout their life.

Any product can be cleaned manually or automatically by washer or disinfectant.

EMS recommends the use of an ISO 15883 compliant automatic washer-disinfector (WD) for optimal effectiveness and service life.

⚠️ EMS handpieces must not be lubricated or otherwise exposed to oil whether directly or indirectly.

AUTOMATIC

CLEAN, DISINFECT AND DRY

! The washer-disinfector must comply with ISO 15883, must have suitable baskets to hold small fragile products and must have rinsing connections with diameter of approximately 16mm for the attachment to product lumen.

Correctly place the product into a suitable rack, connect all lumens to the rinsing connectors and start the automated cleaning.

! Instructions of the manufacturer of the washer disinfectant must also be observed.

⚠️ In case of use of a chemical disinfectant, follow carefully the instructions provided by the disinfection solution manufacturer.

MANUAL

CLEAN

The following process can be used with any EMS product:

Brush the devices in a cleaning solution of 0.5% neodisher¹ in deionized water at 40°C with a suitable soft bristled brush² until all visible residue has been removed.

✓ 15 s For products with lumens, flush all lumens with a spray gun (water jet gun, with a static water pressure of 2 bar) with cold tap water.

✓ 15 min Place the product in 0.5% solution of neodisher¹ with deionized water at 40°C. Make certain that all lumens are filled with cleaning solution (use a syringe if needed).

✓ 15 s Rinse all lumens by flushing with a spray gun (water jet gun, with a static water pressure of 2 bar) with cold deionized water.

✓ 10 s Rinse the whole product under cold running deionized water.

¹ For example: MediClean Dental, Dr. Weigert Hamburg

² Medisafe MED100.33

AUTOMATIC

The following automated process¹ can be used to achieve an A0 Level of 3000:

- ☑ 2 min Pre-clean with cold tap water².
Drain
- ☑ 5 min Clean at 55°C with tap water² and 0.5% cleaning solution.
Drain
- ☑ 3 min Rinse and neutralize with cold demineralized water³.
Drain
- ☑ 2 min Rinse with cold demineralized water³.
Drain
- ☑ 3 min Perform thermal disinfection (final rinse) with deionized water at 93°C minimum.
Drain
- ☑ 20 min Dry at 100°C.
(at least)

MANUAL

 DO NOT use any ultrasound bath cleaning procedure with the handpieces: it may destroy the products.

The following process is intended **ONLY** for EMS PIEZON® Instruments and file holders:

- ☑ 15 s For product with lumens, flush all lumens with a spray gun (water jet gun, with a static water pressure of no more than 2 bar) with cold tap water.
- ☑ 10 min Place products in a metal sieve and sonicate in 0.5% cleaning solution of neodisher⁴ in deionized water at 40°C⁵. Make certain that all lumens are filled with cleaning solution (use a syringe if needed).
- ☑ 15 s Rinse all lumens by flushing with a spray gun (water jet gun, with a static water pressure of no more than 2 bar) with cold deionized water.
- ☑ 10 s Rinse the whole product under cold running deionized water.

DISINFECT

 ASP CIDEX® OPA solution must be used undiluted, within its useful and shelf life, and in compliance with the manufacturer's warnings and instructions for use.

- ☑ 5 min Immerse the product completely into CIDEX® OPA solution at minimum 20°C. Make certain that all lumens are filled with disinfectant solution (use a syringe if needed).
- ☑ 1 min (each rinse) The ASP CIDEX® OPA disinfectant requires a total of three rinses : keep the product totally immersed, and use a large volume⁶ of fresh water.

Do not reuse the water for rinsing or any other purpose. Residues of disinfectant may cause serious side effects.

DRY

Use an air pistol (compressed air) to fully dry the lumen and the whole product until no more residues of water are present (visible or detectable).

INSPECT

 If stains are still visible on the product after cleaning/disinfection, the entire cleaning/disinfection procedure must be repeated. Products with visible damage, chip/flake loss, corrosion or deformation must be disposed of (no further use is allowed). Check also the integrity of O-rings and gaskets and replace if damaged or deformed.

 Verify the product to be fully dry. In case of detection of residues of water, remove these using an air pistol (clean compressed air). Fully dry the lumen and the whole product until no more residues of water are present (visible or detectable).

¹ For example: Miele Professional G 7836 CD with Miele Rack E429

² cold tap water = 16°C +/- 2°C

³ cold demineralized water = 20°C +/- 2°C

⁴ For example: MediClean Dental, Dr. Weigert Hamburg

⁵ For example: Bandelin, Sonorex 1028K, 35kHz

⁶ e.g. 2 gallons

PACK

⚠ Only previously cleaned and disinfected products can be sterilized.

⚠ Effective sterilization can take place only on fully dry products. Ensure each part (internal lumens and any surface) to be perfectly dry before reassembling and packing. Prior to sterilization, the products need to be reassembled for ready use and placed in suitable sterilization packaging.


! AIRFLOW® and PERIOFLOW® Handpieces:

No need for reassembly

! PIEZON® Endochuck:

Reinstall the small O-ring (gasket) first, and then gently screw the retention bolt without tightening.

! PIEZON® Handpieces:

Reinstall the O-ring first, then put the light guide into the nose cap and screw it on the handpiece.

! PIEZON® Instruments with CombiTorque®:

Reinstall the original instrument with its CombiTorque®.

Package your products with a single or double pouch:

- suitable for prevacuum moist heat sterilization,
- compliant with ISO 11607-1 or EN 868,
- resistant to 138°C
- and with adequate steam permeability (e.g. Wipak STERIKING flat rolls Type R43 and R44).

STERILIZE

⚠ Sterilization must be performed immediately after cleaning-disinfection.

The loading pattern defined by the autoclave manufacturer must be followed.

- ⊘ DO NOT exceed**
- the maximum number of sterilization cycles allowed.
 - a sterilization temperature of 138°C and a holding time of 20 min.

⊘ DO NOT use a class N sterilizer to reprocess PIEZON® handpieces.

⊘ DO NOT use hot-air sterilization and radio-sterilization procedures: they destroy the products.

Moist heat sterilization of products shall be performed according to ISO 17665 and in compliance with respective country requirements.

The prevacuum moist heat (steam) process can be used with any EMS product packaged in appropriate single or double pouches:

Parameters for the prevacuum moist heat cycle:

- 3 prevacuum phases
- Pressure of 3 bar¹
- Humidity of 100%
- Temperature of 132°C
- ⌚ 3 min (minimum) Hold time (full cycle)
- ⌚ 20 min (minimum) Dry

⚠ Users must ensure that the reprocessing processes, including resources, materials and personnel, are capable to reach the required results and maintained over time: it is the user's responsibility to ensure that validation of the reprocessing procedures is kept up to date at all times.

STORE

! Store the sterilized products at a temperature of 5°C to 40°C in a:

- dry,
- clean
- and dust-free environment.

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FROM THE VALLEY

► For more than 35 years EMS has developed high-end technologies and protocols in conjunction with the world's most advanced clinicians, enabling the patients to enjoy natural teeth and implants for much longer.

„Primum non nocere“ and minimal invasiveness are not only a must, they are our way of life.

**EMS RESEARCH AND MANUFACTURE CENTER
NYON, SWITZERLAND**

TO THE TOP

MONT BLANC, 4807M

ONLY USE ORIGINAL EMS PRODUCTS

► EMS powders are much less abrasive than traditional mechanical debridement techniques such as rubber cups and handscaling.

We recommend to use AIRFLOW® and PERIOFLOW® first, followed by PIEZON®. GBT: the minimally invasive way for Professional Mechanical Plaque Removal (PMPR).

► As part of the GBT disclosing phase, coloration also minimizes AIRFLOW® abrasion. No more color signifies that all biofilm has been removed: the objective of AIRFLOW® is achieved, to preserve the integrity of dental tissues.

⚠ EMS devices have been designed and tested for exclusive use with EMS prophylaxis powders.

! ONLY USE EMS POWDERS

So-called "EMS compatible" powders in the marketplace may cause damage to oral tissues and the device. See below a damaged Airflow nozzle after use with a "compatible" powder. If EMS products are damaged by using "compatible" powders the product Warranty will be lost.

Damaged handpiece

ONLY USE THE SWISS ORIGINAL EMS INSTRUMENTS!

The EMS Instrument, the Handpiece and the NO PAIN electronic module have to vibrate in harmony like in a trilogy.

Using so-called "compatible" and copy tips from third parties can damage your patients' teeth as well as EMS handpieces.

They were developed and made for each other by EMS.

One PS Instrument is designed for up to 500 treatments, thus reducing the cost of one treatment to only a few cents.

Therefore, there is absolutely no reason to "save" money by using so-called "compatible" tips. Moreover, if the Piezon Handpiece gets damaged (the thread) it will lose its warranty so at the end of the day, "compatible" tips will always cost you more.

HEADQUARTER

NYON, SWITZERLAND

ELECTRO MEDICAL SYSTEMS SA

Chemin de la Vuarpillière 31

CH - 1260 Nyon

www.ems-dental.com/en/contact

www.ems-dental.com

T. +41 22 994 47 00

F. +41 22 994 47 01

WORLDWIDE EMS AFFILIATES

MUNICH, GERMANY

EMS ELECTRO MEDICAL SYSTEMS GMBH

Stahlgruberring 12

DE - 81829 München

T. +49 89 42 71 61 0

F. +49 89 42 71 61 60

E-mail: info@ems-ch.de

FONCINE, FRANCE

EMS ELECTRO MEDICAL SYSTEMS FRANCE SARL

32, route de Pontarlier

FR - 39460 Foncine-le-Haut

T. +33 3 84 51 90 01

F. +33 3 84 51 94 00

E-mail: info@ems-france.fr

MADRID, SPAIN

EMS ELECTRO MEDICAL SYSTEMS ESPAÑA SLU

C/ Tomás Bretón, 50-52 2ª planta

ES - 28045 Madrid

T. +34 91 528 99 89

E-mail: info@ems-espana.com

MILANO, ITALY

EMS ITALIA S.R.L

Via Carlo Pisacane 7B

IT - 20016 Pero (MI)

T. +39 02 3453 8111

E-mail: dental@ems-italia.it

CAMBRIDGESHIRE, ENGLAND

EMS UK LTD SUITE 2

Xenus House - Sandpiper Court

Phoenix Park - Eaton Socon

UK - St Neots, Cambridgeshire

PE19 8EP

T. +44 01480 587260

Email: info@ems-unitedkingdom.com

STOCKHOLM, SWEDEN

EMS NORDIC AB

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SE - 120 30 Stockholm

T. +46 8 899 102

E-mail: info@ems-nordic.se

AMMAN, JORDAN

EMS ELECTRO MEDICAL SYSTEMS JORDAN

No.2 - second floor, 195 Arar (Wadi Saqra) str.

JO - Amman

KATOWICE, POLAND

EMS POLAND SP. Z O.O.

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PL - 40-203 Katowice

T. +48 32 4937060

E-mail: info@ems-poland.com

DALLAS, USA

EMS CORPORATION

11886 Greenville Avenue, Ste 120

US - Dallas, Texas 75243

T. +001 972 690 8382

F. +001 972 690 8981

E-mail: info@ems-na.com

SYDNEY, AUSTRALIA

E.M.S. ELECTRO MEDICAL SYSTEMS OCEANIA

Unit 1, Building 4 - Sydney Corporate Park

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E.M.S. ELECTRO MEDICAL SYSTEMS TRADING

(SHANGHAI) CO., LTD.

24A, Jin Sui Mansion,

N°379 Pudong Nan Rd.

CN - Shanghai 200120

T. +86 21 3363 2323

E-mail: emschina@ems-ch.com

SEOUL, SOUTH KOREA

EMS ELECTRO MEDICAL SYSTEMS KOREA LTD..

13, Beobwon-ro 8-gil, Songpa-gu

T. +82 70 7807 2875

E-mail: emskorea@ems-ch.com

EMS 

WWW.EMS-DENTAL.COM