***MANUALE UTENTE  
USER MANUAL  
ANWENDERHANDBUCH***

# **VITAL TEST**

**Art. 11250**

***ITALIANO***

***ENGLISH***

***DEUTCH***

# INDICE

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# INDEX

## **1. Safety Information**

- 1.1 Warnings
- 1.2 Cautions
- 1.3 Notes
- 1.4 Symbols in the oximeter

## **2. Pulse Oximeter Introduction**

- 2.1 General Introduction
- 2.2 Panels Introduction
  - 2.2.1 Symbols on Screen
  - 2.2.2 Front Panel Buttons
  - 2.2.3 Rear Panel
- 2.3 Connecting Sensor or Cable
- 2.4 Powered by Battery
- 2.5 Accessory List

## **3. Oximeter Operation**

- 3.1 Turn On the Oximeter
- 3.2 Measuring State
  - 3.2.1 Measuring Modes
  - 3.2.2 Trend Graph and Trend Table
  - 3.2.3 Abnormal Measuring State
  - 3.2.4 Data Transfer State
- 3.3 System Menu
  - 3.3.1 System Mode
  - 3.3.2 Patient Type
  - 3.3.3 Alarm Volume
  - 3.3.4 Pulse Volume
  - 3.3.5 Audio Paused (s)
  - 3.3.6 User Maintain
  - 3.3.7 Default Config
  - 3.3.8 Alarm System
  - 3.3.9 SpO2 Alarm setup
  - 3.3.10 PR Alarm setup
  - 3.3.11 Patient ID No. setup
  - 3.3.12 Data Storage
  - 3.3.13 Delete All Data
  - 3.3.14 Exit (Return)

## **4 Alarm**

- 4.1 Alarm Categories and Levels
- 4.2 Alarm conditions
  - 4.2.1 Alarm is off before the first Measurement
  - 4.2.2 Alarm for SpO2 Sensor Unconnected
  - 4.2.3 Alarm for SpO2 Sensor off
  - 4.2.4 Alarm for Low Battery
  - 4.2.5 Higher than Hi Alarm limit
  - 4.2.6 Lower than Lo Alarm limit
  - 4.2.7 Alarm Silence
  - 4.2.8 Turn Off Alarm System
  - 4.2.9 Alarm Priority

## **5 Performance Considerations**

- 5.1 Performance Verification
- 5.2 Oximeter Performance Considerations
- 5.3 Sensor Performance Considerations

## **6 Maintenance**

## **7 Principles of Operation**

- 7.1 Pulse Oximetry Measuring
- 7.2 Functional Versus Fractional Saturation
- 7.3 Measured Versus Calculated Saturation

## **8 Warranty and After-Sales Service**

## **Appendix I**

- A1.1 Classification
- A1.2 Specification
  - A1.2.1 Size and Weight
  - A1.2.2 Environment
  - A1.2.3 Display
  - A1.2.4 Batteries
- A1.3 Parameters

## 1.1 Warnings

Warnings are identified by the WARNING symbol shown above.

Warnings alert the user to potential serious outcomes, such as death, injury, or adverse events to the patient or user.

**WARNING:** Avoid the explosion hazard. Do not use the oximeter in the presence of flammable Anesthetics mixture with air, oxygen, or nitrous oxide.

**WARNING:** Chemicals from a broken LCD display panel are toxic when ingested. Use cautions when the oximeter has a broken display panel.

**WARNING:** Routinely monitor the patient to make sure the oximeter is functioning and the sensor is correctly placed.

**WARNING:** Oximetry measurements and pulse signals can be affected by certain environmental conditions, sensor application errors, and certain patient conditions.

**WARNING:** The use of accessories, sensors, and cables other than those specified may result in increased emission and/or invalid readings of the oximeter.

**WARNING:** Failure to cover the sensor site with opaque material in high ambient light conditions may result in inaccurate measurements.

**WARNING:** Do not silence the audio alarm function, or decrease the audio alarm volume, if patient safety could be compromised.

**WARNING:** The oximeter is a prescription device to be operated only by trained personnel. The oximeter is for attended monitoring only.

**WARNING:** Dispose of batteries in accordance with local ordinances and regulations.

**WARNING:** The oximeter is not defibrillator-proof. However, it may remain attached to the patient throughout defibrillation or while an electrosurgical unit is in use. The measurements may be inaccurate throughout the defibrillation, or use of an electrosurgical unit, and shortly thereafter. To avoid shock, the caregiver should not hold the oximeter while using a defibrillator on a patient.

**WARNING:** Disconnect the oximeter and sensor from the patient throughout magnetic resonance imaging (MRI) scanning. Induced current could potentially cause burns.

**WARNING:** To ensure accurate performance and prevent device failure, do not subject the oximeter to extreme moisture, such as direct exposure to rain. Such exposure may cause inaccurate performance or device failure.

**WARNING:** Do not lift the oximeter by the sensor or extension cable because the cable could disconnect from the oximeter and the oximeter may drop on the patient.

**WARNING:** Do not make any clinical judgment based solely on the oximeter, it is intended only as an adjunct in patient assessment. It must be used in conjunction with clinical signs and symptoms.

**WARNING:** To ensure patient safety, do not place the oximeter in any position that might cause it to fall on the patient.

**WARNING:** As with all medical equipment, carefully route patient cables to reduce the possibility of patient entanglement or strangulation.

**WARNING:** Ensure that the speaker is clear of any obstruction and that the speaker holes are not covered. Failure to do so could result in an inaudible alarm tone.

**WARNING:** Use only manufacturer permitted sensors and extension cables with the oximeter. Other sensors or extension cables may be failed or cause improper monitor performance.

**WARNING:** Pulse Oximeter readings and pulse signals can be affected by certain ambient environmental conditions, sensor application error, and certain patient conditions. See the appropriate sections of the manual for specific safety information.

**WARNING:** Do not let new battery works with battery which will run out; do not use rechargeable batteries works with alkaline batteries.

**WARNING:** Periodically check the battery for corrosion. Take out the batteries from the oximeter if you do not expect to use it within one month.

**WARNING:** If you do not hear the POST pass tone, please do not use the Pulse Oximeter.

**WARNING:** Before using, user should carefully read the applicable user manual of sensor, including warning, cautions and instructions.

**WARNING:** Do not use damaged sensor or extension cables, do not use sensor with exposed optical components.

**WARNING:** Tissue damage may be caused by incorrect application or prolonged measurement duration using the sensor (more than 2 hours). Inspect the sensor periodically according to sensor user manual.

**WARNING:** Do not immerse or wet the sensor, as this may damage the sensor.

**WARNING:** There are no user-serviceable parts inside the oximeter, the cover should only be removed by qualified service personnel.

**WARNING:** Do not spray, pour, spill liquid to the oximeter and its accessories, connector, switch or opening in enclosure, as this may damage the oximeter.

**WARNING:** Before cleaning the oximeter or the sensor, make sure that the equipment is switched off and disconnected from the power line.

## 1.2 Cautions

Cautions are identified by the CAUTION symbol shown above.

Cautions alert the user to exercise care necessary for the safe and effective use of the oximeter.

**CAUTION:** All combinations of equipment must be in compliance with IEC/EN Standard 60601-1-1 systems requirements.

**CAUTION:** The institution should follow local government regulations and recycling instructions regarding disposal or recycling of the batteries and the oximeter components or end of life of the oximeter.

**CAUTION:** Pulse Oximeter will not operate with dead batteries. Install new batteries.

**CAUTION:** The sensor unconnected icon and associated alarm indicate the sensor has disconnected or wire fault. So check the sensor connection and, if necessary, replace the sensor, extension cables or both.

**CAUTION:** During the POST immediately after power-on, confirm all the display segments and icons are shown and the speaker sounds a few seconds tone.

**CAUTION:** When adjusting any menu parameters, the oximeter does not display SpO<sub>2</sub> or PR, but it is still recording.

## 1.3 Notes

**NOTE:**

Notes are identified by the symbol shown above. Notes contain important information that may be overlooked or missed.

**NOTE:**

This device has been tested and found to comply with the requirements of IEC/EN60601-1-2 (International standard for EMC testing of Medical Electrical Equipment, second edition) and the Medical Device Directive 93/42/EEC. These standards are designed to provide reasonable protection against harmful interference in a typical medical installation.

**NOTE:**

Sensor LED light emissions fall within Class 1 level, according to IEC/EN 60825-1:2001. No special safety precautions are required.

**NOTES:**

Normal operation means:

- The oximeter is turned on;
- A sensor is connected to the oximeter;
- The sensor is applied to the patient;
- The patient's SpO<sub>2</sub>, Pulse rate are being reported;
- No error conditions exist.

**NOTE:**

Wash the probe with clean water after disinfecting and sterilizing to remove any remaining solution. The probe can only be reused after being dried thoroughly.

**NOTE:**

Do not disinfect the probe with the means of water boiled.

**NOTE:**

Any residue should be removed from the probe before being disinfected and sterilized, and avoid contacting corrosive solvent. Dip the cable into alcohol or alkalescent solvent for a long time, may reduce the flexibility of the scarfskin of the cable. Also, the connector should not be dipped.


**NOTE:**

After monitoring, disinfect the probe according to the instruction described in the user manual.

**NOTE:**

Manufacturer suggests that the lifetime of the oximeter is 3 years.

## 1.4 Symbols in the oximeter

|                                                                                     |                                                                                                                                                                                      |
|-------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|    | This symbol indicates that the instrument is IEC/EN 60601-1 Type BF equipment.                                                                                                       |
|    | This symbol means 'BE CAREFUL'. Refer to the user manual.                                                                                                                            |
|   | Serial number                                                                                                                                                                        |
|  | CE Mark                                                                                                                                                                              |
|  | Authorised representative in the European community                                                                                                                                  |
|  | Manufacturer                                                                                                                                                                         |
|  | Pursuant to DIRECTIVE 2002/96/EC It indicates that the equipment should be sent to the special agencies according to local regulation for separate collection after its useful life. |

## 2.1 General Introduction

The oximeter is intended for continuous monitoring or spot-checking of functional arterial oxygen saturation (SpO<sub>2</sub>) and pulse rate of single adult, pediatric or neonate patient in hospitals, intra-hospital transport and hospital type facilities.

It displays SpO<sub>2</sub> value, pulse rate value, plethysmogram, bar graph, etc.

The oximeter has been installed SpO<sub>2</sub> module inside. It integrates parameter module function, display and recorder output function. It can be powered by 4 pieces 1.5V AA batteries or 4 pieces 1.2V Ni-H rechargeable AA batteries. It clearly displays all the parameter information on LCD.


Pulse Oximeter

Pulse Oximeter Data Management Software is optional for the oximeter.

## 2.2 Panels Introduction

This section identifies the symbols, controls, displays, and buttons on the front panel of the oximeter.


### 2.2.1 Symbols on Screen

*Buttons and displays*

*Waveform mode*

*Large numeric mode*

Icons on display screen and their meanings:

|                                                                                     |                                                                    |
|-------------------------------------------------------------------------------------|--------------------------------------------------------------------|
| <b>SpO<sub>2</sub></b>                                                              | SpO <sub>2</sub> value display area                                |
| <b>100%</b>                                                                         | Measured SpO <sub>2</sub> %                                        |
| <b>PR</b>                                                                           | Pulse Rate value display area                                      |
| <b>60 bpm</b>                                                                       | Measured Pulse rate (bpm)                                          |
| ↑                                                                                   | Displays when measuring value is higher than the upper alarm limit |
| ↓                                                                                   | Displays when measuring value is lower than the lower alarm limit  |
|  | SpO <sub>2</sub> waveform display                                  |
|  | Pulse amplitude display                                            |
|  | Low battery icon                                                   |
|  | Audio alarm off icon                                               |
|  | Alarm off icon                                                     |
|  | Data storage icon                                                  |
| <b>04:59</b>                                                                        | Time display in Information area: “hour: minute”                   |
| <b>ADU/NEO</b>                                                                      | Patient type in Information area: Adult or Neonate.                |
| <b>ID:99</b>                                                                        | Patient ID in Information area                                     |

|                                                                                   |                                          |
|-----------------------------------------------------------------------------------|------------------------------------------|
|  | SpO <sub>2</sub> sensor unconnected icon |
|  | SpO <sub>2</sub> sensor off              |
|  | Indicates the memory space is full       |
|  | Weak signal icon                         |


#### NOTES:

The icons for sensor unconnected, sensor off or weak signal are displayed in the right part of Information area. Only one of them can be displayed at the same time.


The icon for Patient ID or the icon that indicates the memory space is full are displayed in some part of Information Area. Only one of the two icons can be displayed at the same time.

## 2.2.2 Front Panel Buttons

This section describes the buttons on the front panel of the oximeter. The controls are activated by pressing the button that corresponds to that control. For example, press the Alarm Silence button to control the audio alarm.

*Front Panel buttons*

Turn on or off the oximeter.

On: Press On/Off button and hold it for 1 second.

Off: Press On/Off button and hold it for 2 seconds.

When the oximeter is off, synchronously press On/Off button and Function button for 1 second, the oximeter will enter Data transfer state.

Under the Menu state, press this button to return to measuring state.

During POST, the backlight is not available.

During normal measurement, press this button to turn on or off the backlight.

Alarms that occur throughout the Power-On-Self-Test (POST) can not be silenced.

When Alarm System in menu setup is ON, press Alarm Silence button can turn off audio alarm, the paused period can be set to 30, 60, 90 or 120s. Although the audio alarm is off, the visual alarm is still active. After the paused period is over, the audio alarm is reactivated.

Set the Alarm System to OFF in the menu to turn off the alarm. Pop-up dialog box will display to confirm alarm setting. See details in 3.3.8.

At the menu state, press Up Arrow button to choose different items, and increase the value of some parameters. Press repeatedly to increase a parameter by more than one increment. Press and hold this button for more than 1 second to repeat the increment continuously.

Press this button during measuring state to enter the lately 10 minutes' SpO2 or PR trend graph.

At menu state, press Down Arrow button can choose different items, and decrease the value of some parameters. Press repeatedly to decrease a parameter by more than one decrement. Press and hold the button for more than 1 second to repeat the decrement continuously.

Press this button during measuring state to enter the lately 10 minutes' SpO2 and PR trend table.

During POST, Function button is not available;

Press this button during normal measuring state to enter menu for function choice or setup;


At the menu state, this button is also used as Enter button. Choose one item in menu using cursor button (Up Arrow button and Down Arrow button), and press Function button to confirm, then increase or decrease the value using cursor button.

When the oximeter is off, synchronously press On/Off button and Function button for 1 second, the oximeter will enter Data transfer state.


### **Button Combination**

When the oximeter is off, synchronously press On/Off button and Function button for 1 second, the oximeter will enter Data transfer state.

### 2.2.3 Rear Panel

### 2.3 Connecting Sensor or Cable

*Sensor and Cable Connecting Port*

SpO2 Sensor and cable port is on the top of the oximeter for connecting the SpO2 sensor. A extension cable can be used between the oximeter and the SpO2 sensor. Use only the cable permitted by manufacturer.

Warning symbol

Type BF applied part

Auxiliary output connector

## 2.4 Powered by Battery

The oximeter can be powered by 4 pieces 1.5V LR6 AA alkaline batteries. A new set of batteries can provide up to 48 hours of general operation, or about 24 hours of operation with the backlight and alarm on. The oximeter does not support built-in recharging mode.

The oximeter also can be powered by 4 pieces 1.2V Ni-H rechargeable batteries.

### Battery Installation

To install batteries:

1. Make sure the oximeter has been turned off.
2. Pull the battery compartment cover downward toward the bottom of the oximeter, and remove it.
3. Install 4 AA batteries according to the indications on battery compartment of the oximeter.
4. Replace the battery compartment cover.

### Low Battery Icon

The Low battery icon displays and alarm is given when few minutes operation remains available. After few minutes operation, the oximeter will turn off automatically. Replace the batteries.

Low Battery Icon

## 2.5 Accessory List

Standard configuration including:


| QUANTITY | DESCRIPTION                |
|----------|----------------------------|
| 1        | Pulseoximeter              |
| 4        | Battery AA 1,5 V (IEC LR6) |
| 1        | SpO <sub>2</sub> sensor    |
| 1        | User manual                |
| 1        | Bag                        |

### 3.1 Turn On the Oximeter

Power-On-Self-Test (POST) verifies the oximeter works properly and is safe to use. Proper operation of the oximeter is verified every time it is turned on as described in the procedure. The verification procedure POST needs approximately few seconds to complete. When the oximeter is turned on, POST automatically tests the circuitry and functions of the oximeter.

Press the **On/Off** button for 1s to turn on the oximeter.

- ◆ At first the LOGO is shown.

- ◆ If the POST is successfully finished, the oximeter sounds a tone, then enter the main interface.

If there is error during POST, the following error codes will display on the screen:

| Error code | Indication                                 |
|------------|--------------------------------------------|
| Error 01   | Indicates error for Low battery            |
| Error 02   | Indicates error for SpO <sub>2</sub> board |
| Error 03   | Indicates error for Main control board     |

## 3.2 Measuring State

### 3.2.1 Measuring Modes

There two measuring modes which are Waveform mode and Large numeric mode. The default configuration is Waveform mode.

#### Waveform Mode


At normal measuring state, oximeter can measure arterial oxygen saturation and pulse rate, display oxygen saturation level and symbol (%SpO<sub>2</sub>) and PR on interface. Besides, it can also display pulse bar graph and Plethysmogram.

*Waveform mode*

#### Large Numeric Mode

The oximeter can display SpO<sub>2</sub>, oxygen saturation unit (%), PR, pulse rate unit (bpm) in large numeric mode.


*Large numeric mode*

### 3.2.2 Trend Graph and Trend Table

During normal measuring state, press Up Arrow button to enter lately SpO2 or PR trend graph, and press Down Arrow button to enter lately 10 minutes' SpO2 and PR trend table. Shift the pages by pressing Up Arrow or Down Arrow button.

Trend graph:

Display SpO2 and PR trend graph

Trend table:

| TREND TABLE |      |     |
|-------------|------|-----|
| TIME        | SpO2 | PR  |
| 20:00:06    | 100  | 66  |
| 20:00:00    | 99   | 68  |
| 19:59:54    | ---  | --- |
| 19:59:48    | ---  | --- |
| 19:59:42    | 98   | 62  |

Display SpO2 and PR trend table

### 3.2.3 Abnormal Measuring State

If the SpO2 sensor does not connect to the oximeter, it will give Medium alarm, and display  in the Information area.

If the SpO2 sensor falls off from the finger, it will give Medium alarm, and display  in the Information area.

At menu state or trend state, without operation for 30 seconds, the oximeter will return to measuring state.

In measuring state, if there is not measuring data and without operation for 10 minutes, the oximeter will turn off automatically.

At Data transfer state, if the oximeter does not receive responsible signals for 10 minutes, it will turn off automatically.

### 3.2.4 Data Transfer State

Set Data Storage in menu to ON, the measuring value will be stored in the oximeter. The SpO2 and PR information can be transferred from oximeter to Pulse Oximeter Data Management Software.


Procedure for data transferring:

- After the measurement and storage are all finished, turn off the oximeter;
- Connect the cable to oximeter and computer for communicating with the Pulse Oximeter Data Management Software;
- Synchronously press On/Off button and Function button, after POST, the oximeter enters Data Transfer State automatically.

## 3.3 System Menu

Press Function button to see the following main menu of oximeter, select items by Up/Down button, and confirm it by Function button.

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#### Alarm Setup:

|                                 |            |
|---------------------------------|------------|
| <b>Alarm System</b>             | <b>ON</b>  |
| <b>SpO<sub>2</sub> Hi Alarm</b> | <b>100</b> |
| <b>SpO<sub>2</sub> Lo Alarm</b> | <b>90</b>  |
| <b>PR Hi Alarm</b>              | <b>120</b> |

|                    |           |
|--------------------|-----------|
| <b>PR Lo Alarm</b> | <b>50</b> |
| <b>Return</b>      |           |

#### Storage Setup:

|                        |            |
|------------------------|------------|
| <b>Patient ID No.</b>  | <b>1</b>   |
| <b>Data Storage</b>    | <b>OFF</b> |
| <b>Delete All Data</b> |            |
| <b>Return</b>          |            |

#### NOTE:

The SpO<sub>2</sub> Hi Alarm and SpO<sub>2</sub> Lo Alarm stand for the upper and lower alarm limit of SpO<sub>2</sub>; The PR Hi Alarm and PR Lo Alarm stand for the upper and lower alarm limit of PR. If user changes the default value for Lo Alarm or Hi Alarm, after restarting the oximeter, the value will resume to the default value for corresponding patient type.

There are some items in the menu, hereinafter all the items will be introduced.

### 3.3.1 System Mode

There are two items for selecting:

Waveform mode 

Large numeric mode 

Then confirm the selection by Function button.

### 3.3.2 Patient Type

The Patient Type can be set to different measuring mode: Adu for Adult or Pediatric patient, or Neo for Neonate patient.

Set Patient Type to Adu or Neo, and confirm it by Function button.

### 3.3.3 Alarm Volume

The Alarm Volume is used to adjust alarm volume, which ranges from 1 to 5.

When the Alarm System is ON, if Low alarm, Medium alarm or High alarm occurs, the oximeter sounds beep.

### 3.3.4 Pulse Volume

User can turn on or off the pulse volume by Pulse Volume, and change volume level to 1, 2, 3, 4, 5 or OFF. Press Function button to enter setup state, then use Up Arrow or Down Arrow button to choose, then confirm it by Function button.

The oximeter implements variable pulse tone whose frequency varies according to the saturation.


### 3.3.5 Audio Paused (s)

Set the paused period for audio alarm to 30, 60, 90 or 120 seconds.

When the Alarm System is ON, press Alarm Silence button can turn off the audio alarm, the paused period is set by the Audio Paused (s).

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Enter User Maintain menu through password (819):


**Enter Password**

8  9

**Confirm      Cancel**


If the password is wrong, the following dialog box will pop up:

**Wrong Password!!**  
**Please Retry.**

**EXIT**

If the password is right, the following menu will display:


**Language** **ENG**

**Time Setup >>**

**Return**

- Language: user can select language to display.
- Time Setup >>: select this item, the following interface displays:


**2008 -  - 01**

**18 : 56 : 59**

**Confirm      Cancel**

### 3.3.7 Default Config

Choose this item to resume factory default configuration. A dialog box pops up:


Factory Default Configuration is as following list:

|                   |                                                                                     |
|-------------------|-------------------------------------------------------------------------------------|
| System Mode:      |  |
| Patient Type:     | ADU                                                                                 |
| Alarm System:     | ON                                                                                  |
| Alarm Volume:     | 3                                                                                   |
| Pulse Volume:     | 3                                                                                   |
| Audio Paused (s): | 60                                                                                  |
| SpO2 Hi Alarm:    | 100                                                                                 |
| SpO2 Lo Alarm:    | 90                                                                                  |
| PR Hi Alarm:      | 120                                                                                 |
| PR Lo Alarm:      | 50                                                                                  |
| Patient ID No.:   | 1                                                                                   |
| Data Storage:     | OFF                                                                                 |

### 3.3.8 Alarm System

Set the Alarm System to ON or OFF to turn on or off the alarm system.

If the Alarm system is set to OFF, a dialog box pops up as follows:

If the **Alarm System** is **ON**, when alarm occurs, the oximeter gives visual alarm and audio alarm.

Press **Alarm Silence** button can pause the alarm system for seconds (the paused period can be set to 30, 60, 90 or 120s by user, see 3.3.5), the Audio alarm off icon displays. But the visual alarm is still active. For example, if the measured  $\text{SpO}_2$  value is higher than **SpO<sub>2</sub> Hi Alarm** or lower than **SpO<sub>2</sub> Lo Alarm**, there will be  $\uparrow$  or  $\downarrow$  icon display on screen, and the  $\text{SpO}_2$  or PR character will flash.

If the **Alarm system** is set to **OFF**, all the audio alarm and visual alarm are turned off.

#### **WARNING**

When the Alarm system is off, the oximeter will not give alarm prompt. In order to avoid endangering the patient's life, the user should use this function cautiously.

### 3.3.9 $\text{SpO}_2$ Alarm setup

User can choose **SpO<sub>2</sub> Hi Alarm** and **SpO<sub>2</sub> Lo Alarm** in menu to adjust  $\text{SpO}_2$  alarm limit. Press **Up Arrow** button or **Down Arrow** button to increase or decrease alarm limit.

In default configuration setup, the **SpO<sub>2</sub> Hi Alarm** and **SpO<sub>2</sub> Lo Alarm** for **Neo** mode are 95 and 90; while they are 100 and 90 for **Adu** mode.

Set the  $\text{SpO}_2$  alarm limits as following procedure:

- ◆ Choose **SpO<sub>2</sub> Hi Alarm** in the menu, press **Function** button to enter setup. The **SpO<sub>2</sub> Hi Alarm** box will change from real line box to broken line box. The adjustable range for upper limit of  $\text{SpO}_2$  is from 1 plus the lower limit of  $\text{SpO}_2$  to 100.
- ◆ Press **Up Arrow** or **Down Arrow** button to increase or decrease values.
- ◆ Choose **SpO<sub>2</sub> Lo Alarm** in the menu, press **Function** button to set up. The **SpO<sub>2</sub> Lo Alarm** box will change from real line box to broken line box. The adjustable

range for the lower limit of SpO<sub>2</sub> Alarm is from 0 to the upper limit of SpO<sub>2</sub> Alarm minus 1.

- ◆ Press **Up Arrow** or **Down Arrow** button to increase or decrease values.
- ◆ **SpO<sub>2</sub> Hi Alarm** is always higher than **SpO<sub>2</sub> Lo Alarm** at least 1%.
- ◆ Press **Function** button, confirm the alarm range setup.
- ◆ Press **On/Off** button to exit menu, and return to measuring state.

### 3.3.10 PR Alarm setup

User can use **PR Hi Alarm** and **PR Lo Alarm** in menu to adjust pulse rate alarm limits.

In default configuration setup, the **PR Hi Alarm** and **PR Lo Alarm** for **Neo** mode are 200 and 100; while they are 120 and 50 for **Adu** mode.

Set the PR limits as following procedure:


- ◆ Choose **PR Hi Alarm** in the menu, press **Function** button to enter setup. The **PR Hi Alarm** box changes from real line to broken line. The adjustable range of the upper limit of PR Alarm is from 1 plus the lower limit of PR Alarm to 250.
- ◆ Press **Up Arrow** or **Down Arrow** button to increase or decrease values.
- ◆ Choose **PR Lo Alarm** in menu, press **Function** button enter setup. The **PR Lo Alarm** box changes from real line to broken line. The adjustable range for the lower limit of PR Alarm is from 0 to the upper limit of PR Alarm minus 1.
- ◆ Press **Function** button, confirm the alarm range setup.
- ◆ **Hi Alarm** is always higher than **Lo Alarm** at least 1 bpm.
- ◆ Press **On/Off** button to exit menu, and return to measuring state.

### 3.3.11 Patient ID No. setup

The oximeter can support 100 patient ID, and 300 hours data storage time.

When entering menu, press **Function** button to set ID (valid range is from 1 to 100). ID display box on interface will change from real line to broken line.

After choosing ID, press **Function** button to confirm the setup. If the ID exists, dialog box will pop up for confirming to cover the former data.

### 3.3.12 Data Storage

Choose **Data Storage**, and set to **ON**, then the measuring data can be stored.

During data storage, patient ID can not be changed. If the user wants to change ID, he should change **Data Storage** to **OFF**, then setup new ID.


Data stored in the oximeter can be exported through HPC100 Pulse Oximeter Data Management Software. Please refer to 3.2.4 for Data transfer procedure.

When the memory space is full, a icon  displays in Information Area. Meanwhile **Data Storage** changes to **OFF** automatically. Restart the oximeter, a dialog box pops up, user should confirm to delete all the data.

### 3.3.13 Delete All Data

**Delete All Data** is for deleting all the stored data. Choose this item by **Function** button, a dialog box pops up as follows:

If you choose **YES** for deleting all the data, the deleting progress shows:

### 3.3.14 Exit (Return)

Exit menu by **Exit** in menu.

Return to the previous menu by **Return** in the menu.

## **4 Alarm**

### **4.1 Alarm Categories and Levels**

#### **Alarm Categories:**

The oximeter's alarms can be classified into two categories: physiological alarms and technical alarms.

1. Physiological alarms

Physiological alarms, also called patient status alarms, are triggered by a monitored parameter value that violates set alarm limits or an abnormal patient condition.

2. Technical alarms

Technical alarms, also called system status alarms, are triggered by a device malfunction or a patient data distortion due to improper operation or system problems.

#### **Alarm Levels**

By severity, the oximeter's alarms levels can be classified into three categories: high level alarms, medium level alarms and low level alarms.

1. High level alarms

Indicate that the patient is in a life threatening situation and an emergency treatment is demanded.

2. Medium level alarms

The patient's vital signs appear abnormal or the oximeter system status appears abnormal, indicate that prompt operator response is required.

3. Low level alarms

The patient's vital signs appear abnormal or the oximeter system status appears abnormal, indicate that operator awareness is required.

The levels for both technical alarms and physiological alarms are predefined and can not be changed by user.

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- Character flash
- Alarm tone

High level alarms: character flashes quickly and sounds triple + double + triple +double beep;

Medium level alarms: character flashes slowly and sounds triple beep;

Low level alarms: character display constantly and sounds a single beep.

## 4.2 Alarm conditions

### 4.2.1 Alarm is off before the first Measurement

Before the first measurement, the alarm system is configured to be off. At this time, if the SpO<sub>2</sub> sensor is unconnected or the sensor is off, the oximeter will not give alarm.

### 4.2.2 Alarm for SpO<sub>2</sub> Sensor Unconnected

When the SpO<sub>2</sub> sensor is disconnected, the oximeter gives Medium alarm. The icon  displays in Information area.

It displays --- in SpO<sub>2</sub>, PR value area of LCD, and gives Medium alarm. (Make sure the **Alarm System** in menu is ON.)

### 4.2.3 Alarm for SpO<sub>2</sub> Sensor off

When the SpO<sub>2</sub> sensor falls off from the finger, the oximeter will give Medium alarm, and the icon  displays in Information Area.

It displays --- in SpO<sub>2</sub>, PR value area of LCD, and gives Medium alarm. (Make sure the **Alarm System** in menu is ON.)

### 4.2.4 Alarm for Low Battery

When battery is too low, the oximeter gives Medium alarm for Low battery.

After the Low battery alarm occurs, the oximeter can still be operated for few minutes before it turns off automatically.

The Low battery icon  displays on LCD, and gives Medium alarm. (Make sure the **Alarm System** in menu is ON.)

## 4.2.5 Higher than Hi Alarm limit

If the measured SpO<sub>2</sub> or PR value is higher than the **Hi Alarm** (upper alarm limit), the oximeter gives High alarm.

Here take PR for example:

When the measured PR value is higher than the setup **PR Hi Alarm**, the oximeter gives High alarm (Make sure the **Alarm System** in menu is **ON**). A ↑ icon displays near PR, which indicates that the measured value is higher than that of **PR Hi Alarm**, it will synchronously flash with PR value.

## 4.2.6 Lower than Lo Alarm limit

If the measured SpO<sub>2</sub> or PR value is lower than the **Lo Alarm** (lower alarm limit), the oximeter gives High alarm.

Here take SpO<sub>2</sub> for example:

When the measured SpO<sub>2</sub> value is lower than the setup **SpO<sub>2</sub> Lo Alarm**, the oximeter gives Low SpO<sub>2</sub> alarm. (Make sure the **Alarm System** in menu is **ON**.)

A ↓ icon displays near SpO<sub>2</sub>, which indicates the measured value is lower than that of **SpO<sub>2</sub> Lo Alarm**, it will synchronously flash with SpO<sub>2</sub> value.

## 4.2.7 Alarm Silence

If the **Alarm System** in menu is **ON**, press **Alarm Silence** button, the audio alarm will be off for the paused period set by user, but the visual alarm still active.

When the audio alarm is off, press **Alarm Silence** button to reactivate the audio alarm function.

## 4.2.8 Turn Off Alarm System

After the **Alarm system** is turned off, the oximeter can not give visual and audio alarms except for Low battery icon alarm.

Set the **Alarm system** to **ON**, the alarm system will be active. It will give audio alarm and visual alarm if alarm occurs.

## 4.2.9 Alarm Priority

Only one kind of alarm can be given at once. For example, if Medium alarm and High alarm occur at the same time, High alarm will take priority.

If the pulse beep and audio alarm sound at the same time, the oximeter will only give alarm sound.

## 5.1 Performance Verification

Qualified service personnel are responsible for performance verification procedures before the oximeter is used for the first time in a clinical setting.

## 5.2 Oximeter Performance Considerations

There are some patient conditions that can affect the oximeter's measurements.

- ◆ **Dysfunctional Hemoglobins**

Dysfunctional hemoglobins, such as carboxyhemoglobin, methemoglobin, and sulfhemoglobin, are unable to carry oxygen. SpO<sub>2</sub> readings may appear normal; however, a patient may be hypoxic because less hemoglobin is available to carry oxygen. Further assessment beyond pulse oximeter is recommended.

- ◆ **Anemia**

Anemia causes decreased arterial oxygen content. Although SpO<sub>2</sub> readings may appear normal, an anemic patient may be hypoxic. Correcting anemia can improve arterial oxygen content. The oximeter may fail to provide SpO<sub>2</sub> if hemoglobin levels fall below 5 gm/dl.

- ◆ **Saturation**

The oximeter displays saturation level between 1% and 100%.

- ◆ **Pulse rate**

The oximeter displays pulse rate between 30 and 254 beats per minute (bpm). The sensor accuracy ranges do not apply to pulse rates above 254 bpm.

## 5.3 Sensor Performance Considerations

Inaccurate measurements can be caused by:

- ◆ Incorrect application of the sensor.
- ◆ Placement of the sensor on an extremity with a blood pressure cuff, arterial catheter, or intravascular line.
- ◆ Excessive patient activity.
- ◆ Intravascular dyes, such as indocyanine green or methylene blue.
- ◆ Externally applied coloring, such as nail polish or pigmented cream.
- ◆ Failure to cover the sensor site with opaque materials in high ambient light conditions.
- ◆ Venous pulsation.
- ◆ Dysfunctional hemoglobin.
- ◆ Low perfusion.

Loss-of-pulse signal occurs for the following reasons:

- The sensor is applied too tightly.
- Defibrillation.
- A blood pressure cuff is inflated on the same extremity as the one with the sensor attached.
- There is arterial occlusion proximal to the sensor.
- Poor peripheral perfusion.
- Loss of pulse/cardiac arrest.

To use the sensor:

- Select an appropriate sensor.
- Apply the sensor as directed, and observe all warnings and cautions presented in the sensor user manual.
- Clean and remove any substances, such as nail polish, from the application site.
- Periodically check to ensure that the sensor remains properly positioned on the patient.

High ambient light sources that can interfere with the performance of the sensor are:

- Surgical lights (especially those with a xenon light source).
- Bilirubin lamps.
- Fluorescent lights.
- Infrared heating lamps.
- Direct sunlight.

To prevent interference from ambient light, ensure that the sensor is properly applied, and cover the sensor site with opaque material.

If interference due to patient activity presents a problem, try one or more of the following to correct the problem:

- Verify that the sensor is properly and securely applied.
- Move the sensor to another site.
- Use an adhesive to the sensor.
- Use a new the sensor with fresh adhesive backing.
- Keep the patient still, if possible.

**WARNING**

Before cleaning the oximeter or the sensor, make sure that the oximeter is switched off. The oximeter does not require calibration.

If service is necessary, contact qualified service personnel or your local manufacturer representative.

Before using the oximeter, do the following:

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If you find any damage on the oximeter, stop using the oximeter on patient, and contact the biomedical engineer of the hospital or Customer service immediately.

#### Periodic Safety Checks

It is recommended that the following checks be performed every 24 months:

- ◆ Inspect the devices for mechanical and functional damage
- ◆ Inspect the relevant labels for legibility

All the checks that need to open the oximeter should be performed by qualified customer service technician. The safety and maintenance check can be conducted by persons from this company. You can obtain the material about the customer service contract from the local company's office.

If the hospital or agency that is responding to using the oximeter does not follow a satisfactory maintenance schedule, the oximeter may become invalid, and the human health may be endangered.

#### Cleaning

You can surface-clean and disinfect the oximeter and sensor.

To surface-clean the oximeter:


- ◆ Use a soft cloth dampened with either a commercial, nonabrasive cleaner, or a solution of 70% alcohol in water.
- ◆ Lightly wipe the surfaces of the oximeter.

To disinfect the oximeter:

- ◆ Use a soft cloth saturated with a solution of 10% chlorine bleach in tap water.

Pulse Oximeter adopts non-invasive double wavelength to measure  $\text{SpO}_2$  and PR. It can perform spot measuring and continuous measuring for a short time.

The system is consisted of Central Processing Unit, Signal Collection, Signal Input, Data Output, Display and User Input module, shown as follows:

Communication between the oximeter and external devices is through RS-232 interface.

## 7.1 Pulse Oximetry Measuring

The oximeter uses oximetry to measure functional oxygen saturation in the blood. Pulse oximetry works by applying sensor to a pulsating arteriolar vascular bed, such as a finger or toe. The sensor contains a dual light source and a photonic detector.

Bone, tissue, pigmentation, and venous vessels normally absorb a constant amount of light over time. The arteriolar bed normally pulsates and absorbs variable amounts of light during the pulsations. The ratio of light absorbed is translated into a measurement of functional oxygen saturation ( $\text{SpO}_2$ ). Because a measurement of  $\text{SpO}_2$  is dependent upon light from the sensor, excessive ambient light can interfere with this measurement. Pulse oximetry is based on two principles:

- ◆ Oxyhemoglobin and deoxyhemoglobin differ in their absorption of red and infrared light (spectrophotometry).
- ◆ The volume of arterial blood in tissue (hence light absorption by the blood) changes during the pulse (plethysmography).

The oximeter determines  $\text{SpO}_2$  by passing red and infrared light into an arteriolar bed and measuring changes in light absorption during the pulsatile cycle. Red and infrared low-

voltage light-emitting diodes (LED) serve as light sources; a photonic diode serves as the photo detector.

Because oxyhemoglobin and deoxyhemoglobin differ in light absorption, the amount of red and infrared light absorbed by blood is related to hemoglobin oxygen saturation. To identify the oxygen saturation of arterial hemoglobin, the oximeter uses the pulsatile nature of arterial flow.

During systole, a new pulse of arterial blood enters the vascular bed, and blood volume and light absorption increase. During diastole, blood volume and light absorption reach their lowest point.

The oximeter bases its SpO<sub>2</sub> measurements on the difference between maximum and minimum absorption (measurements at systole and diastole). By doing so, it focuses on light absorption by pulsatile arterial blood, eliminating the effects of non-pulsatile absorbers such as tissue, bone and venous blood.

### **Wavelength**

The sensor contains LEDs that emit red light at a wavelength of approximately 660 nm and infrared light at a wavelength of approximately 900 nm.

## **7.2 Functional Versus Fractional Saturation**

This oximeter measures functional saturation-oxygenated hemoglobin expressed as a percentage of the hemoglobin that can transport oxygen. It does not detect significant amounts of dysfunctional hemoglobin, such as carboxyhemoglobin or methemoglobin.

In contrast, hemoximeter such as the IL482 report fractional saturation-oxygenated hemoglobin expressed as a percentage of all measured hemoglobin, including measured dysfunctional hemoglobins.

To compare functional saturation measurements to those from an instrument that measures fractional saturation, fractional measurements must be converted as follows:

$$\text{Functional saturation} = \frac{\text{Fractional saturation}}{\times 100 - (\% \text{Carboxyhemoglobin} + \% \text{Methemoglobin})}$$

## **7.3 Measured Versus Calculated Saturation**

When saturation is calculated from a blood gas partial pressure of oxygen (PO<sub>2</sub>), the calculated value may differ from the SpO<sub>2</sub> measurement of a pulse oximeter. This usually occurs because the calculated saturation was not appropriately corrected for the effects of variables that shift the relationship between PO<sub>2</sub> and pH, the partial pressure of carbon dioxide (PCO<sub>2</sub>), 2,3-DPG, and fetal hemoglobin.

1) DURATION: the product is covered by a 12-month (twelve) guarantee.

2) EXCLUSION: the guarantee excludes:

- a) faults originating from non-performance of scheduled maintenance and due to negligence or improper use of the product by the User;
- b) periodic checks and maintenance;
- c) repair or replacement of parts subject to wear that are fragile or have an unpredictable life, unless they prove faulty at the time of their delivery;
- d) faults caused by personnel or originating from technical personnel transfer and transportation;
- e) faults or damage originating from improper use or operating errors;
- f) faults or damage originating from contamination in the water and air infeed systems, extraordinary chemical agents or electrical events;
- g) faults or damage originating from the use of detergents, disinfecting agents, sterilizing fluids or processes that are not described explicitly in the instruction manual for operation and maintenance;
- h) natural color change of plastic components.

3) INSTALLATION AND TESTING REPORT: an essential prerequisite for guarantee validity is the return of the installation and testing report after all its sections have been filled in by the User and the Fitter. The report must be returned within 15 days after installation, otherwise the guarantee becomes null and void.

4) LIMITS: the guarantee entitles its holder to repairs or replacement of faulty components free of charge. The right to replace the entire machine is excluded. As for components either applied or added to TECNO-GAZ products or by third parties, having their own guarantee certificate, the conditions, limits and exclusions indicated in the relevant certificates apply. Under no circumstances shall carrying out of one or more than one repair during the guarantee term change the guarantee expiry date.

5) CLAIMS: in the event of the Purchaser making a claim relating to the application of the guarantee and the quality or conditions of the delivered equipment, the said Purchaser shall in no way suspend and/or delay payments.

6) DISPUTES: if a dispute arises on the application and interpretation of this Certificate of guarantee, the Court of PARMA (Italy) shall be the court of jurisdiction, regardless of the location in which the purchase agreement was stipulated.

7) DEROGATIONS: any special derogations from the guarantee terms described in this section shall not imply the granting to the Purchaser of any right and shall be considered as granted for a specific case only.

8) MISCELLANEOUS: for any other issue, which is not covered by this Certificate of Guarantee, the rules described in the Civil Code and existing relevant Laws in the Republic of Italy apply.

## Appendix I Specification

### A1.1 Classification

|                      |                                        |
|----------------------|----------------------------------------|
| Type of Protection   | Internally powered equipment           |
| EMC Compliance       | Class B                                |
| Degree of Protection | Type BF-Applied part                   |
| Ingress Protection   | IPX2                                   |
| Mode of operation    | Continuous measuring and spot checking |

### A1.2 Specification

#### A1.2.1 Size and Weight

|        |                              |
|--------|------------------------------|
| Size   | 160 (L)×70 (W)×37.6 (H) (mm) |
| Weight | 165 (g) (without battery)    |

#### A1.2.2 Environment

|             |                |
|-------------|----------------|
| Temperature |                |
| Working     | 5 °C ~ 40 °C   |
| Storage     | -20 °C ~ 55 °C |

|          |                          |
|----------|--------------------------|
| Humidity |                          |
| Working  | 25% ~ 80% (No coagulate) |
| Storage  | 25% ~ 93% (No coagulate) |

|                       |                    |
|-----------------------|--------------------|
| Atmospheric pressure  |                    |
| Working               | 860 hPa ~ 1060 hPa |
| Transport and Storage | 700 hPa ~ 1060 hPa |

#### A1.2.3 Display

|                  |                                                              |
|------------------|--------------------------------------------------------------|
| Screen Type      | 128×64 dot-matrix LCD, with white LED backlight              |
| Big Numeric Mode | SpO <sub>2</sub> , PR and Bar graph displayed                |
| Waveform Mode    | SpO <sub>2</sub> , PR, Bar graph and Plethysmogram displayed |

#### A1.2.4 Batteries

|                           |                                          |
|---------------------------|------------------------------------------|
| Power supply              | 4 Batteries                              |
| Battery life is typically | 48 hours                                 |
| Standard match            | 4 pieces 1.5V AA alkaline batteries;     |
| Compatible with           | 4 pieces 1.2V Ni-H rechargeable battery; |
|                           | Can not recharge with battery build-in   |

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**TECNO-GAZ S.p.A.**